

- **New low- cost entry level instrument . Developed for users with basic requirements at an affordable price.**



The versatile all-round moisture determination solution

- Sturdy construction, perfect for all laboratory and production environments
- Weighing technology which meets all international metrological regulations
- Weighing range 52 g / readability of balance 1 mg
- High temperature range of 170°C, temperature increments 1°C, results displayed to .01%
- Practical construction to enable easy cleaning
- High contrast, bright vacuum fluorescent display
- One method memory storage only
- RS232 interface
- Easy-clean design
- Halogen or infra-red heater options

Simple, fast and accurate

- Variety of start-up and end-point control selections for optimum control of sample treatment
- Language independent prompts operated by symbol keys
- Easy access sample holder

Quality Assurance

- Results displayed in GLP guideline format. Numerical and statistical analysis available
- Real-time results shown in display
- Download latest software through the internet

Food industry

- Special dark infra-red version without glass is available for the food industry